

PTPLAD1 Polyclonal Antibody

Catalog No: #29868



Package Size: #29868-1 50ul #29868-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTPLAD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PTPLAD1 (NP_057479.2).
Other Names	HACD3; B-IND1; BIND1; HSPC121; PTPLAD1; 3-hydroxyacyl-CoA dehydratase 3
Accession No.	Swiss-Prot#:Q9P035NCBI Gene ID:51495
Uniprot	Q9P035
GeneID	51495;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only