

PCF11 Polyclonal Antibody

Catalog No: #29869



Package Size: #29869-1 50ul #29869-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PCF11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PCF11 (NP_056969.2).
Other Names	PCF11; pre-mRNA cleavage complex 2 protein Pcf11
Accession No.	Swiss-Prot#:O94913NCBI Gene ID:51585
Uniprot	O94913
GeneID	51585;
Calculated MW	180kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PCF11 at 1:1000 dilution.

Background

The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor IIm. The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage.

Note: This product is for in vitro research use only