

ZCCHC10 Polyclonal Antibody

Catalog No: #29870



Package Size: #29870-1 50ul #29870-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZCCHC10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ZCCHC10 (NP_001287745.1).
Other Names	ZCCHC10; zinc finger CCHC-type containing 10
Accession No.	Swiss-Prot#:Q8TBK6NCBI Gene ID:54819
Uniprot	Q8TBK6
GeneID	54819;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only