

OGFOD1 Polyclonal Antibody

Catalog No: #29873

Package Size: #29873-1 50ul #29873-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OGFOD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human OGFOD1 (NP_060703.3).
Other Names	OGFOD1; TPA1; prolyl 3-hydroxylase OGFOD1
Accession No.	Swiss-Prot#:Q8N543NCBI Gene ID:55239
Uniprot	Q8N543
GeneID	55239;
Calculated MW	63kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of HeLa cells, using OGFOD1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only