

LIN7C Polyclonal Antibody

Catalog No: #29874



Package Size: #29874-1 50ul #29874-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LIN7C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LIN7C (NP_060832.1).
Other Names	LIN7C; LIN-7-C; LIN-7C; MALS-3; MALS3; VELI3; protein lin-7 homolog C
Accession No.	Swiss-Prot#:Q9NUP9NCBI Gene ID:55327
Uniprot	Q9NUP9
GeneID	55327;
Calculated MW	25kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LIN7C at 1:1000 dilution.

Background

Note: This product is for in vitro research use only