

# CPXM1 Polyclonal Antibody

Catalog No: #29878

Package Size: #29878-1 50ul #29878-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CPXM1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human CPXM1
Other Names	CPXM1; CPX1; CPXM; probable carboxypeptidase X1
Accession No.	Swiss-Prot#:Q96SM3NCBI Gene ID:56265
Uniprot	Q96SM3
GeneID	56265;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB□1:500 - 1:2000

## Background

This gene likely encodes a member of the carboxypeptidase family of proteins. Cloning of a comparable locus in mouse indicates that the encoded protein contains a discoidin domain and a carboxypeptidase domain, but the protein appears to lack residues necessary for carboxypeptidase activity.

Note: This product is for in vitro research use only