

GPR108 Polyclonal Antibody

Catalog No: #29879



Package Size: #29879-1 50ul #29879-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GPR108 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GPR108 (NP_001073921.1).
Other Names	GPR108; LUSTR2; protein GPR108
Accession No.	Swiss-Prot#:Q9NPR9NCBI Gene ID:56927
Uniprot	Q9NPR9
GeneID	56927;
Calculated MW	80kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GPR108 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only