

ZNF471 Polyclonal Antibody

Catalog No: #29882

Package Size: #29882-1 50ul #29882-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF471 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF471 (NP_065864.2).
Other Names	ZNF471; ERP1; Z1971; zinc finger protein 471
Accession No.	Swiss-Prot#:Q9BX82NCBI Gene ID:57573
Uniprot	Q9BX82
GeneID	57573;
Calculated MW	90kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF471 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only