

usp29 Polyclonal Antibody

Catalog No: #29885



Package Size: #29885-1 50ul #29885-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	usp29 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human usp29 (NP_065954.1).
Other Names	USP29; HOM-TES-84/86; ubiquitin specific peptidase 29
Accession No.	Swiss-Prot#:Q9HBJ7NCBI Gene ID:57663
Uniprot	Q9HBJ7
GeneID	57663;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Background

Note: This product is for in vitro research use only