

FAM65A Polyclonal Antibody

Catalog No: #29887



Package Size: #29887-1 50ul #29887-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM65A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human FAM65A (NP_001180451.1).
Other Names	RIPOR1; FAM65A; protein FAM65A
Accession No.	Swiss-Prot#:Q6ZS17NCBI Gene ID:79567
Uniprot	Q6ZS17
GeneID	79567;
Calculated MW	140kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse kidney, using FAM65A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only