

VCPIP1 Polyclonal Antibody

Catalog No: #29889



Package Size: #29889-1 50ul #29889-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	VCPIP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human VCPIP1 (NP_079330.2).
Other Names	VCPIP1; DUBA3; VCIP135; deubiquitinating protein VCIP135
Accession No.	Swiss-Prot#:Q96JH7NCBI Gene ID:80124
Uniprot	Q96JH7
GeneID	80124;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using VCPIP1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only