

TNKS2 Polyclonal Antibody

Catalog No: #29890



Package Size: #29890-1 50ul #29890-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TNKS2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TNKS2 (NP_079511.1).
Other Names	TNKS2; ARTD6; PARP-5b; PARP-5c; PARP5B; PARP5C; TANK2; TNKL; pART6; tankyrase-2
Accession No.	Swiss-Prot#:Q9H2K2NCBI Gene ID:80351
Uniprot	Q9H2K2
GeneID	80351;
Calculated MW	126kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TNKS2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only