

STARD3NL Polyclonal Antibody

Catalog No: #29891

Package Size: #29891-1 50ul #29891-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STARD3NL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human STARD3NL (NP_114405.1).
Other Names	STARD3NL;MENTHO
Accession No.	Uniprot:O95772GeneID:83930
Uniprot	O95772
GeneID	83930
Calculated MW	30kDa
SDS-PAGE MW	26KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded Human colon carcinoma using STARD3NL Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat brain using STARD3NL Rabbit pAb.



Western blot analysis of extracts of various cell lines, using STARD3NL antibody.

Background

This gene encodes a late-endosomal protein that contains a conserved MENTAL (MLN64 N-terminal) domain. The encoded protein binds cholesterol molecules and may play a role in endosomal cholesterol transport through interactions with metastatic lymph node protein 64 (MLN64).

Note: This product is for in vitro research use only