

DHX57 Polyclonal Antibody

Catalog No: #29895



Package Size: #29895-1 50ul #29895-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DHX57 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human DHX57 (NP_945314.1).
Other Names	DHX57; DDX57; DExH-box helicase 57
Accession No.	Swiss-Prot#:Q6P158NCBI Gene ID:90957
Uniprot	Q6P158
GeneID	90957;
Calculated MW	150-155kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DHX57 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only