

TDRD9 Polyclonal Antibody

Catalog No: #29896

Package Size: #29896-1 50ul #29896-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TDRD9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TDRD9 (NP_694591.2).
Other Names	TDRD9; C14orf75; HIG-1; NET54; tudor domain containing 9
Accession No.	Swiss-Prot#:Q8NDG6NCBI Gene ID:122402
Uniprot	Q8NDG6
GeneID	122402;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Background

Note: This product is for in vitro research use only