

SYNPR Polyclonal Antibody

Catalog No: #29899



Package Size: #29899-1 50ul #29899-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SYNPR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SYNPR (NP_653243.1).
Other Names	SYNPR; SPO; synaptoporin
Accession No.	Swiss-Prot#:Q8TBG9NCBI Gene ID:132204
Uniprot	Q8TBG9
GeneID	132204;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only