

NOTUM Polyclonal Antibody

Catalog No: #29903

Package Size: #29903-1 50ul #29903-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOTUM Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human NOTUM (NP_848588.3).
Other Names	NOTUM; hNOTUM; palmitoleoyl-protein carboxylesterase NOTUM
Accession No.	Swiss-Prot#:Q6P988NCBI Gene ID:147111
Uniprot	Q6P988
GeneID	147111;
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NOTUM at 1:1000 dilution.

Background

Note: This product is for in vitro research use only