

SPC24 Polyclonal Antibody

Catalog No: #29904



Package Size: #29904-1 50ul #29904-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SPC24 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SPC24 (NP_872319.1).
Other Names	SPC24; SPBC24; kinetochore protein Spc24
Accession No.	Swiss-Prot#:Q8NBT2NCBI Gene ID:147841
Uniprot	Q8NBT2
GeneID	147841;
Calculated MW	22kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SPC24 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only