

RBM46 Polyclonal Antibody

Catalog No: #29908



Package Size: #29908-1 50ul #29908-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RBM46 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RBM46 (NP_659416.1).
Other Names	RBM46;CT68
Accession No.	Uniprot:Q8TBY0GeneID:166863
Uniprot	Q8TBY0
GeneID	166863
Calculated MW	56kDa
SDS-PAGE MW	65KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RBM46 antibody.

Note: This product is for in vitro research use only