

CFAP53 Polyclonal Antibody

Catalog No: #29909

Package Size: #29909-1 50ul #29909-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CFAP53 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CFAP53 (NP_659457.2).
Other Names	CFAP53;CCDC11;HTX6
Accession No.	Uniprot:Q96M91GeneID:220136
Uniprot	Q96M91
GeneID	220136
Calculated MW	62kDa
SDS-PAGE MW	62kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of HeLa cells using CFAP53 Rabbit pAb.



Western blot analysis of extracts of various cell lines, using CFAP53 antibody.

Background

This gene belongs to the CFAP53 family. It was found to be differentially expressed by the ciliated cells of frog epidermis and in skin fibroblasts from human. Mutations in this gene are associated with visceral heterotaxy-6, which implicates this gene in determination of left-right asymmetric patterning.

Note: This product is for in vitro research use only