

RNF152 Polyclonal Antibody

Catalog No: #29911



Package Size: #29911-1 50ul #29911-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	RNF152 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	A synthetic peptide of human RNF152 (NP_775828.1).
Other Names	RNF152; ring finger protein 152
Accession No.	Swiss-Prot#:Q8N8N0NCBI Gene ID:220441
Uniprot	Q8N8N0
GeneID	220441;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Background

Note: This product is for in vitro research use only