

LNK2 Polyclonal Antibody

Catalog No: #29913



Package Size: #29913-1 50ul #29913-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LNK2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LNK2 (NP_699202.1).
Other Names	LNK2; PDZRN1; ligand of Numb protein X 2
Accession No.	Swiss-Prot#:Q8N448NCBI Gene ID:222484
Uniprot	Q8N448
GeneID	222484;
Calculated MW	76kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LNK2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only