

ZNF763 Polyclonal Antibody

Catalog No: #29916



Package Size: #29916-1 50ul #29916-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF763 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZNF763 (NP_001012771.1).
Other Names	ZNF763; ZNF; ZNF440L; zinc finger protein 763
Accession No.	Swiss-Prot#:Q0D2J5NCBI Gene ID:284390
Uniprot	Q0D2J5
GeneID	284390;
Calculated MW	46kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF763 at 1:1000 dilution.



Western blot analysis of extracts of various cell lines, using ZNF763 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only