

ATP11C Polyclonal Antibody

Catalog No: #29917



Package Size: #29917-1 50ul #29917-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ATP11C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ATP11C (NP_775965.2).
Other Names	ATP11C; ATPIG; ATPIQ; phospholipid-transporting ATPase IG
Accession No.	Swiss-Prot#:Q8NB49NCBI Gene ID:286410
Uniprot	Q8NB49
GeneID	286410;
Calculated MW	157kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ATP11C at 1:1000 dilution.

Background

Note: This product is for in vitro research use only