

HRNR Polyclonal Antibody

Catalog No: #29920



Package Size: #29920-1 50ul #29920-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HRNR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human HRNR (NP_001009931.1).
Other Names	HRNR; FLG3; S100A16; S100a18; hornerin
Accession No.	Swiss-Prot#:Q86YZ3NCBI Gene ID:388697
Uniprot	Q86YZ3
GeneID	388697;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF 1:50 - 1:200

Images



Immunofluorescence analysis of HeLa cells using HRNR at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only