

ZSCAN1 Polyclonal Antibody

Catalog No: #29922



Package Size: #29922-1 50ul #29922-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZSCAN1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZSCAN1 (NP_872378.3).
Other Names	ZSCAN1; MZF-1; ZNF915; zinc finger and SCAN domain containing 1
Accession No.	Swiss-Prot#:Q8NBB4NCBI Gene ID:284312
Uniprot	Q8NBB4
GeneID	284312;
Calculated MW	45kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of HepG2 cells, using ZSCAN1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only