

TMEM100 Polyclonal Antibody

Catalog No: #29927



Package Size: #29927-1 50ul #29927-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TMEM100 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic synthetic peptide of human TMEM100 (NP_001093110.1).
Other Names	TMEM100; transmembrane protein 100
Accession No.	Swiss-Prot#:Q9NV29NCBI Gene ID:55273
Uniprot	Q9NV29
GeneID	55273;
Calculated MW	14kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TMEM100 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only