

LSM14A Polyclonal Antibody

Catalog No: #29933

Package Size: #29933-1 50ul #29933-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LSM14A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human LSM14A
Other Names	LSM14A; C19orf13; FAM61A; RAP55; RAP55A; protein LSM14 homolog A
Accession No.	Swiss-Prot#:Q8ND56NCBI Gene ID:26065
Uniprot	Q8ND56
GeneID	26065;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of L929 cells using LSM14A at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only