

CALN1 Polyclonal Antibody

Catalog No: #29934

Package Size: #29934-1 50ul #29934-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CALN1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human CALN1 (NP_001017440.1).
Other Names	CALN1; CABP8; calneuron 1
Accession No.	Swiss-Prot#:Q9BXU9NCBI Gene ID:83698
Uniprot	Q9BXU9
GeneID	83698;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

This gene encodes a protein with high similarity to the calcium-binding proteins of the calmodulin family. The encoded protein contains two EF-hand domains and potential calcium-binding sites. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only