

## NLRC5 Polyclonal Antibody

Catalog No: #29943

Package Size: #29943-1 50ul #29943-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NLRC5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human NLRC5
Other Names	NLRC5;CLR16.1;NOD27;NOD4
Accession No.	Uniprot:Q86WI3GeneID:84166
Uniprot	Q86WI3
GeneID	84166
Calculated MW	202kDa
SDS-PAGE MW	202kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of mouse lung, using NLRC5 antibody.

## Background

This gene encodes a member of the caspase recruitment domain-containing NLR family. This gene plays a role in cytokine response and antiviral immunity through its inhibition of NF-kappa-B activation and negative regulation of type I interferon signaling pathways.

---

Note: This product is for in vitro research use only