

# IQSEC1 Polyclonal Antibody

Catalog No: #29945



Package Size: #29945-1 50ul #29945-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                          |
|-----------------------|--------------------------------------------------------------------------|
| Product Name          | IQSEC1 Polyclonal Antibody                                               |
| Host Species          | Rabbit                                                                   |
| Clonality             | Polyclonal                                                               |
| Isotype               | IgG                                                                      |
| Purification          | Affinity purification                                                    |
| Applications          | WB                                                                       |
| Species Reactivity    | Mouse,Rat                                                                |
| Immunogen Description | Recombinant fusion protein of human IQSEC1 (NP_055684.3).                |
| Other Names           | IQSEC1; ARF-GEP100; ARFGEP100; BRAG2; GEP100; IQ motif and Sec7 domain 1 |
| Accession No.         | Swiss-Prot#:Q6DN90NCBI Gene ID:9922                                      |
| Uniprot               | Q6DN90                                                                   |
| GeneID                | 9922;                                                                    |
| Calculated MW         | 108kDa                                                                   |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.        |
| Storage               | Store at -20°C                                                           |

## Application Details

WB 1:500 - 1:1000

## Images



Western blot analysis of extracts of various cell lines, using IQSEC1 at 1:1000 dilution.

## Background

---

Note: This product is for in vitro research use only