

TFAP4 Polyclonal Antibody

Catalog No: #29946

Package Size: #29946-1 50ul #29946-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TFAP4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Rat
Immunogen Description	Synthetic Peptide of human TFAP4
Other Names	TFAP4; AP-4; bHLHc41; transcription factor AP-4
Accession No.	Swiss-Prot#:Q01664NCBI Gene ID:7023
Uniprot	Q01664
GeneID	7023;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:200 - 1:2000IF□1:50 - 1:200


Images




Immunofluorescence analysis of C6 cells using TFAP4 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using TFAP4 at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts from normal (control) and TFAP4 knockout (KO) HeLa cells, using TFAP4 at 1:500 dilution.



Western blot analysis of extracts of various cell lines, using TFAP4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only