

PSD3 Polyclonal Antibody

Catalog No: #29955

Package Size: #29955-1 50ul #29955-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PSD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human PSD3.
Other Names	PSD3; EFA6D; EFA6R; HCA67; PH and SEC7 domain-containing protein 3
Accession No.	Swiss-Prot#:Q9NYI0NCBI Gene ID:23362
Uniprot	Q9NYI0
GeneID	23362;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse brain, using PSD3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only