

# SFSWAP Polyclonal Antibody

Catalog No: #29959



Package Size: #29959-1 50ul #29959-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SFSWAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SFSWAP (NP_004583.2).
Other Names	SFSWAP; SFRS8; SWAP; splicing factor SWAP homolog
Accession No.	Swiss-Prot#:Q12872NCBI Gene ID:6433
Uniprot	Q12872
GeneID	6433;
Calculated MW	150kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using SFSWAP at 1:1000 dilution.

## Background

This gene encodes a human homolog of Drosophila splicing regulatory protein. This gene autoregulates its expression by control of splicing of its first two introns. In addition, it also regulates the splicing of fibronectin and CD45 genes. Two transcript variants encoding different isoforms have been identified.

---

Note: This product is for in vitro research use only