

## ZNF711 Polyclonal Antibody

Catalog No: #29960

Package Size: #29960-1 50ul #29960-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF711 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZNF711.
Other Names	ZNF711; CMPX1; MRX97; ZNF4; ZNF5; ZNF6; Zfp711; dJ75N13.1; zinc finger protein 711
Accession No.	Swiss-Prot#:Q9Y462NCBI Gene ID:7552
Uniprot	Q9Y462
GeneID	7552;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF711 at 1:1000 dilution.

## Background

This gene encodes a zinc finger protein of unknown function. It bears similarity to a zinc finger protein which acts as a transcriptional activator. This gene lies in a region of the X chromosome which has been associated with mental retardation.

---

Note: This product is for in vitro research use only