

TNKS Polyclonal Antibody

Catalog No: #29962



Package Size: #29962-1 50ul #29962-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TNKS Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human TNKS (NP_003738.2).
Other Names	TNKS; ARTD5; PARP-5a; PARP5A; PARPL; TIN1; TINF1; TNKS1; pART5; tankyrase
Accession No.	Swiss-Prot#:O95271NCBI Gene ID:8658
Uniprot	O95271
GeneID	8658;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TNKS at 1:1000 dilution.

Background

Note: This product is for in vitro research use only