

GTF2B Polyclonal Antibody

Catalog No: #29969



Package Size: #29969-1 50ul #29969-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GTF2B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human GTF2B (NP_001505.1).
Other Names	GTF2B;TF2B;TFIIB
Accession No.	Uniprot:Q00403GeneID:2959
Uniprot	Q00403
GeneID	2959
Calculated MW	35kDa
SDS-PAGE MW	32KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using GTF2B antibody.




Immunohistochemistry of paraffin-embedded human lung cancer using GTF2B Rabbit pAb.




Immunohistochemistry of paraffin-embedded human esophageal using GTF2B Rabbit pAb.



Immunofluorescence analysis of L929 cells using GTF2B Polyclonal Antibody.



Immunohistochemistry of paraffin-embedded mouse testis using GTF2B Rabbit pAb.



Immunohistochemistry of paraffin-embedded human colon using GTF2B Rabbit pAb.

Background

This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.

Note: This product is for in vitro research use only