

RPP40 Polyclonal Antibody

Catalog No: #29973



Package Size: #29973-1 50ul #29973-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RPP40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RPP40 (NP_001273061).
Other Names	RPP40; RNASEP1; bA428J1.3; ribonuclease P/MRP subunit p40
Accession No.	Swiss-Prot#:O75818NCBI Gene ID:10799
Uniprot	O75818
GeneID	10799;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RPP40 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only