

PMF1 Polyclonal Antibody

Catalog No: #29976



Package Size: #29976-1 50ul #29976-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PMF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human PMF1 (NP_009152.2).
Other Names	PMF1; polyamine-modulated factor 1
Accession No.	Swiss-Prot#:Q6P1K2NCBI Gene ID:11243
Uniprot	Q6P1K2
GeneID	100527963;11243;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PMF1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only