

ARFIP2 Polyclonal Antibody

Catalog No: #29979

Package Size: #29979-1 50ul #29979-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARFIP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ARFIP2 (NP_036534.1).
Other Names	ARFIP2; POR1; arfaptin-2
Accession No.	Swiss-Prot#:P53365NCBI Gene ID:23647
Uniprot	P53365
GeneID	23647;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded human breast cancer using ARFIP2 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using ARFIP2 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only