

## PRPF6 Polyclonal Antibody

Catalog No: #29980

Package Size: #29980-1 50ul #29980-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PRPF6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human PRPF6 (NP_036601.2).
Other Names	PRPF6; ANT-1; ANT1; C20orf14; Prp6; RP60; SNRNP102; TOM; U5-102K; hPrp6; pre-mRNA-processing factor 6
Accession No.	Swiss-Prot#:O94906NCBI Gene ID:24148
Uniprot	O94906
GeneID	24148;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of HT-29 cells, using PRPF6 at 1:1000 dilution.

## Background

The protein encoded by this gene appears to be involved in pre-mRNA splicing, possibly acting as a bridging factor between U5 and U4/U6 snRNPs in formation of the spliceosome. The encoded protein also can bind androgen receptor, providing a link between transcriptional activation and splicing.

---

Note: This product is for in vitro research use only