

PRPF6 Polyclonal Antibody

Catalog No: #29980

Package Size: #29980-1 50ul #29980-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRPF6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human PRPF6 (NP_036601.2).
Other Names	PRPF6; ANT-1; ANT1; C20orf14; Prp6; RP60; SNRNP102; TOM; U5-102K; hPrp6; pre-mRNA-processing factor 6
Accession No.	Swiss-Prot#:O94906NCBI Gene ID:24148
Uniprot	O94906
GeneID	24148;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Images



Western blot analysis of extracts of HT-29 cells, using PRPF6 at 1:1000 dilution.

Background

The protein encoded by this gene appears to be involved in pre-mRNA splicing, possibly acting as a bridging factor between U5 and U4/U6 snRNPs in formation of the spliceosome. The encoded protein also can bind androgen receptor, providing a link between transcriptional activation and splicing.

Note: This product is for in vitro research use only