

GAPVD1 Polyclonal Antibody

Catalog No: #29981

Package Size: #29981-1 50ul #29981-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GAPVD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human GAPVD1 (NP_056450.2).
Other Names	GAPVD1; GAPEX5; GAPex-5; RAP6; GTPase activating protein and VPS9 domains 1
Accession No.	Swiss-Prot#:Q14C86NCBI Gene ID:26130
Uniprot	Q14C86
GeneID	26130;
Calculated MW	165kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:100

Images



Immunofluorescence analysis of L929 cells using GAPVD1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using GAPVD1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only