

## CNNM4 Polyclonal Antibody

Catalog No: #29982

Package Size: #29982-1 50ul #29982-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CNNM4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human CNNM4.
Other Names	CNNM4;ACDP4
Accession No.	Uniprot:Q6P4Q7GeneID:26504
Uniprot	Q6P4Q7
GeneID	26504
Calculated MW	Refer to figures
SDS-PAGE MW	87kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB□1:500 - 1:2000IF□1:50 - 1:100

## Images



Western blot analysis of extracts of various cell lines, using CNNM4 antibody.

## Background

This gene encodes a member of the ancient conserved domain containing protein family. Members of this protein family contain a cyclin box motif and have structural similarity to the cyclins. The encoded protein may play a role in metal ion transport. Mutations in this gene are associated with Jalili syndrome which consists of cone-rod dystrophy and amelogenesis imperfecta.

---

Note: This product is for in vitro research use only