

FCF1 Polyclonal Antibody

Catalog No: #29983



Package Size: #29983-1 50ul #29983-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FCF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human FCF1.
Other Names	FCF1;Bka;C14orf111;CGI-35;UTP24
Accession No.	Uniprot:Q9Y324GeneID:51077
Uniprot	Q9Y324
GeneID	51077
Calculated MW	Refer to figures
SDS-PAGE MW	23kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of 293T cells, using FCF1 antibody.



Immunohistochemistry of paraffin-embedded Human esophageal using FCF1 Rabbit pAb.



Immunofluorescence analysis of U-2 OS cells using FCF1 antibody.



Immunofluorescence analysis of C6 cells using FCF1 antibody.

Background

Essential protein involved in pre-rRNA processing and 40S ribosomal subunit assembly.

Note: This product is for in vitro research use only