

ARL6IP4 Polyclonal Antibody

Catalog No: #29984

Package Size: #29984-1 50ul #29984-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARL6IP4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ARL6IP4 (NP_057722.3).
Other Names	ARL6IP4; SFRS20; SR-25; SRp25; SRp37; ADP ribosylation factor like GTPase 6 interacting protein 4
Accession No.	Swiss-Prot#:Q66PJ3NCBI Gene ID:51329
Uniprot	Q66PJ3
GeneID	51329;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded human lung using ARL6IP4 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast cancer using ARL6IP4 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only