

AJAP1 Polyclonal Antibody

Catalog No: #29991

Package Size: #29991-1 50ul #29991-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AJAP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human AJAP1.
Other Names	AJAP1; MOT8; SHREW-1; SHREW1; adherens junction-associated protein 1
Accession No.	Swiss-Prot#:Q9UKB5NCBI Gene ID:55966
Uniprot	Q9UKB5
GeneID	55966;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using AJAP1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only