

OSTC Polyclonal Antibody

Catalog No: #29993



Package Size: #29993-1 50ul #29993-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	OSTC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic synthetic peptide of human OSTC (NP_067050.1).
Other Names	OSTC; DC2; oligosaccharyltransferase complex subunit OSTC
Accession No.	Swiss-Prot#:Q9NRP0NCBI Gene ID:58505
Uniprot	Q9NRP0
GeneID	58505;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using OSTC at 1:1000 dilution.

Background

Note: This product is for in vitro research use only