

UBE2O Polyclonal Antibody

Catalog No: #29994

Package Size: #29994-1 50ul #29994-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBE2O Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human UBE2O (NP_071349.3).
Other Names	UBE2O; E2-230K; ubiquitin conjugating enzyme E2 O
Accession No.	Swiss-Prot#:Q9C0C9NCBI Gene ID:63893
Uniprot	Q9C0C9
GeneID	63893;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC 1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded human colon carcinoma using UBE2O at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast cancer using UBE2O at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using UBE2O at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only