

# ATOH8 Polyclonal Antibody

Catalog No: #29999



Package Size: #29999-1 50ul #29999-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | ATOH8 Polyclonal Antibody                                         |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB                                                                |
| Species Reactivity    | Mouse,Rat                                                         |
| Immunogen Description | A synthetic peptide of human ATOH8.                               |
| Other Names           | ATOH8; HATH6; bHLHa21; protein atonal homolog 8                   |
| Accession No.         | Swiss-Prot#:Q96SQ7NCBI Gene ID:84913                              |
| Uniprot               | Q96SQ7                                                            |
| GeneID                | 84913;                                                            |
| Calculated MW         | 35kDa                                                             |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ATOH8 at 1:1000 dilution.

## Background

---

Note: This product is for in vitro research use only